

ELSODENT

Clinical and more biocompatible products

Healbond

Light-curing bonding agent



Description

- Healbond is an advanced strong adhesive system. The surface to be treated has to be etched beforehand.

Features

- Healbond has good wetting properties, to be used on slightly moist dentin surfaces (wet-bonding-technique).
- Contains benzalkonium chloride to eliminate residual bacteria (no post-operative sensitivities).
- Due to its hydrophilic properties, HEALBOND deeply penetrates into the tubulis, to provide a strong retention and an almost complete elimination of residual bacteria.
- Releases fluoride.
- Tensile strength after 10 minutes: ~ 17 MPa.

Indications

- To bond light-curing composites in dental cavities.
- To desensitize cervical regions and free dentin areas.

Reference & Presentation

B-5	5 ml bottle.
-----	--------------

MD: class IIA. Certifying body: SNCH. Manufacturer: G-Pharma



Phone: +33 (0)1 30 37 75 75 - elsodent@elsodent.com www.elsodent.com
ELSODENT, 35 Avenue du Gros Chêne - BP 50072 - HERBLAY - 95612 CERGY PONTOISE CEDEX - FRANCE