

Notre Dam

Light-curing liquid dam

TEGDMA, HEMA, BISGMA and BPA-free



Features

- Light-curing resin whose thixotropy allows an easy setting of the material.
- Low exothermic reaction to avoid harmful reaction of the mucosa.
- Versatile in use.
- Blue for a good contrast.

Indications

- **Bleaching**
 - Protection of soft tissues during bleaching process.
- **Endodontic**
 - To insure a waterproof joint of the rubber dam around a tooth during an endodontic session.
 - To built temporary tooth walls during endodontic sessions, to hold the irrigating liquid in place.
- **Prosthodontic**
 - To fill up undercuts before impression, to avoid tearing off the material.
- **Restaurative dentistry**
 - To hold Teflon® film in place when luting aesthetic veneers, or making a class II composite involving the contact point.



1 - Before reconstruction of dentinal wall.



2 - Dentinal walls built.



Used as waterproof joint around a tooth during an endodontic session.



Gum protection during whitening.

References & Presentations

ND- 4	2 x 2g syringes + 10 black applicator tips. Blue shade.
LTL	100 black needle tips. Gauge 20.

MD: class IIA. Certifying body: SNCH. Manufacturer: G-Pharma

